



Dr.-Ing. RITTER Apparatebau GmbH & Co KG
Coloniastr. 19-23
D-44892 Bochum
Germany

Tel.: +49-(0)234 - 92293-0
Fax: +49-(0)234 - 92293-50
Email: mailbox@ritter.de
Website: www.ritter.de

Divisions:
Plastics engineering & manufacturing
Instruments engineering
Measuring instruments

Content

1. General Information	4
1.1. This Documentation	4
1.2. Positioning the sensor	4
1.3. Treatment of Gases prior to Analysis.....	4
1.4. Additional tube connection for flushing of the casing	5
1.5. Casing Types	6
1.6. Dimensions	6
1.7. Assignment of Sensor Types to possible Sensor Combinations	7
1.8. Calibrations	14
1.9. Measurement of Humid Gases	15
2. Disadvantages of electrochemical (EC) sensors in comparison to the optical gas sensors »RITTER MultiGas«	16
3. Scope of Supply	16
4. Setup	17
5. Infrared Sensor.....	18
5.1. Description	18
5.2. Applications.....	18
5.3. Specifications	18
5.4. Table of Standard Measurement Ranges	20
5.5. Table of Detection Limits (= 3 σ) in Percent of Full Scale (FS)	21
5.6. Recalibrations	21
6. Ultraviolet Sensors.....	22
6.1. Description	22
6.2. Applications.....	22
6.3. Specifications NDUV Sensor	22
6.4. Specifications H ₂ S Sensor.....	24
6.5. Specifications UVRAS Sensor.....	26
6.6. Table of Standard Measurement Ranges	28
6.7. Table of Detection Limits (= 3 σ) in Percent of Full Scale	28
6.8. Recalibrations	29
7. Options (installed inside of the sensor casing)	30
7.1. Oxygen Sensor	30
7.2. Pressure Sensor	32
Pressure Sensor for IR Sensors	32
Pressure Sensor for UV Sensors.....	32
7.3. Humidity Sensor	33

7.4.	Automatic Zero-point Adjustment by Air	35
7.5.	Automatic Zero-point Adjustment by Nitrogen	36
7.6.	LCD-Display for Sensor Casing.....	37
7.7.	Analog Voltage Output	37
7.8.	Heating and Thermostating of Sensor Casing	38
7.9.	Tube Port for Flushing of Casing	39
7.10.	Filter	39
8.	Display for Sensor Casing Type 1.....	41
8.1.	Stand Bracket / Carrying Handle	41
8.2.	Startup Phase	41
8.3.	Measurement Display Screen.....	42
8.4.	Measurement Display Refresh	42
8.5.	Default Settings.....	42
8.6.	Menu Navigation – Description of Menu Items	42
8.7.	Menu Navigation – Overview of Menu Items.....	45

1. General Information

1.1. This Documentation

This documentation has been prepared with the greatest possible care. However, it is intended **for specialist personnel** who are familiar with the principles of gas measurement technology. Basic questions about measurements of gas concentrations, cross-sensitivities etc. can therefore not be answered with this documentation. If necessary, please consult the relevant special literature and/or qualified personnel.

1.2. Positioning the sensor

The sensor must not be used:

- in potentially explosive atmospheres
- in flow direction after a MilliGascounter

Additional recommendations:

- The gas to be analysed should have a stable flow rate between 5 and 300 ltr/h.
- If an external pump is used it should have no pressure fluctuations and no back pressure

1.3. Treatment of Gases prior to Analysis

- The gas must be **free of particles** in order to prevent a blockage of the sensor cuvette. RITTER offers an optional inline particle filter (see par. 7.10 "Filter"). However, it's the user's responsibility to check whether this filter is suitable for the intended application. Furthermore, the filter must be checked regularly und replaced if necessary.
- **Condensation must not occur** inside of the sensor. Condensation can be prevented by:
 - Drying of the gas to be analysed to control the dew point to a value below that of the sensor's ambient temperature. This can be done by using an electric cooler.
 - Condensation trap in the gas line in front of the sensor
 - Heating the sensor and/or sensor cuvettes
(see par. 7.8 "Heating and Thermostating of Sensor Casing")
- The gas **must not contain aggressive components** which could damage the sensor materials unless it has been expressly designed for this gas(es). The used materials are described ...
 - ... in par. 5.3 for NDIR sensors and
 - ... in par. 6.3 ff for NDUV sensors.

1.4. Additional tube connection for flushing of the casing

The gas lines inside the casing and the measuring cuvette are gas-tight by means of O-rings and other gas-tight connections. However, as with all connections, a leakage rate, however slight, cannot be ruled out. With an additional tube connection for flushing of the casing, an accumulation of the measuring gas inside the casing can be prevented.


In the case of **toxic or aggressive gases**, a suction line can be connected via this tube connection, which creates an under-pressure inside the housing thus preventing the gas from escaping from the casing.

In case of **ignitable gases (methane, hydrogen, etc.)**, a pressure or suction line can be connected via this tube connection creating an over- or under-pressure inside the casing thus preventing the accumulation of an ignitable gas mixture.


Please note: H₂S can be perceived as an unpleasant smell even in the ppb range. Even with an O-ring seal, micro-leaks can occur, which can be extracted by flushing the housing.

1.5. Casing Types

Type 1



Type 2



1.6. Dimensions

	Type 1	Type 2
A	171 mm	444 mm
B	290 mm ¹⁾	305 mm
C	86 mm	145 mm
Weight, approx. ²⁾	2 ⁺ kg	Empty: 7.5 kg Incl. 9 x IR: 13 kg
Connections gas in-/outlet	PVDF screw-type tube connection for tube Ø 4 _i / 6 _o mm	
Interfaces	Standard: USB Options: RS232, CANbus, CANopen	
Power supply	24 V DC, incl. plug-in power supply unit Input: 100-240 V-AC, 50/60 Hz	

¹⁾ May be larger with cuvettes for ppm measurement ranges

²⁾ Depending on type and number of built-in sensor modules

1.7. Assignment of Sensor Types to possible Sensor Combinations

Part 1: Pre-configured sensors, incl. casing type 1 with gas connection fittings, data interface, plug-in power supply - ready for use

Part 2: Pre-configured sensors, incl. casing type 2 with gas connection glands, data interface, fittings, data interface, plug-in power supply - ready for use

Part 3: Individually combinable sensor modules, composition according to customer specification, delivery incl. housing in suitable size according to number and type of modules - ready for use

Part 4: Casings for sensor modules, assembled according to customer's specification

Part 5: Built-in options mounted inside of the casing

Part 6: Accessories / options

Part 1: Pre-configured Sensors Incl. casing type 1 with gas connection fittings, data port, power supply Ready for use					
Sensor Type: RMS- xxx	Article No.	Group of Gases IR	Number of Detectable Gases in this Group	Group of Gases UV	Number of Detectable Gases in this Group
xxx = mono IR1	2678	CO ₂ CO N ₂ O CH ₄ C _n H _m * CF ₄ SF ₆	1		
xxx = duo IR2	2742	CO ₂ CO N ₂ O CH ₄ C _n H _m * CF ₄ SF ₆	2		

* Calibration with Propane

Sensor Type: RMS- xxx	Article No.	Group of Gases IR	Number of Detectable Gases in this Group	Group of Gases UV	Number of Detectable Gases in this Group
xxx = trio IR3	2743	CO ₂ CO N ₂ O CH ₄ C _n H _m * CF ₄ SF ₆	3		
xxx = mono UV1	2749			SO ₂ NO ₂ O ₃ Cl ₂ ≤0.5%	1
xxx = mono UV1 resist	2763			SO ₂ Cl ₂ ≤30%	1
xxx = duo UV2	2766			SO ₂ ≤0.5% NO ₂ O ₃ Cl ₂ ≤0.5%	2
xxx = duo IR1+UV1	2797	CO ₂ CO N ₂ O CH ₄ C _n H _m * CF ₄ SF ₆	1	SO ₂ NO ₂ O ₃ Cl ₂ ≤5%	1
xxx = trio IR1+UV2		CO ₂ CO N ₂ O CH ₄ C _n H _m * CF ₄ SF ₆	1	SO ₂ NO ₂ O ₃ Cl ₂ ≤0.5%	2

* Calibration with Propane

Part 2: Pre-configured Sensors Incl. casing type 2 with gas connection fittings, data port, power supply Ready for use					
Sensor Type: RMS- xxx	Article No.	Group of Gases IR	Number of Detectable Gases in this Group	Group of Gases UV	Number of Detectable Gases in this Group
xxx = mono UV1 H₂S ≤ 5.000 ppm	2672			H ₂ S ≤ 5.000 ppm	1
xxx = mono UV1 H₂S ≤ 1%	2855			H ₂ S ≤ 1%	1
xxx = UVRAS	2812			SO ₂ NO ₂ NO	3
xxx = duo IR1 + H₂S ≤ 5000 ppm	2959	CO ₂ CO N ₂ O CH ₄ C _n H _m * CF ₄ SF ₆	1	H ₂ S ≤5000ppm	1
xxx = duo IR1 + [H₂S ≤1% / NO]	2960	CO ₂ CO N ₂ O CH ₄ C _n H _m * CF ₄ SF ₆	1	H ₂ S ≤1% NO	1

* Calibration with Propane

Part 3: Individually combinable sensor modules Composition according Customer Specification Delivery incl. casing of suitable size according to number and type of modules Ready for use

Sensor Type: RMS- xxx	Article No.	Group of Gases IR	Number of Detectable Gases in this Group	Group of Gases UV	Number of Detectable Gases in this Group
xxx = Mod IR1	2813	CO ₂ CO N ₂ O CH ₄ C _n H _m * CF ₄ SF ₆	1		
xxx = Mod IR2	2814	CO ₂ CO N ₂ O CH ₄ C _n H _m * CF ₄ SF ₆	2		
xxx = Mod IR3	2815	CO ₂ CO N ₂ O CH ₄ C _n H _m * CF ₄ SF ₆	3		

* Calibration with Propane

Sensor Type: RMS- xxx	Article No.	Group of Gases IR	Number of Detectable Gases in this Group	Group of Gases UV	Number of Detectable Gases in this Group
xxx = Mod UV1	2830			SO ₂ NO ₂ O ₃ Cl ₂ ≤0.5%	1
xxx = Mod UV1 resist				SO ₂ Cl ₂ ≤30	1
xxx = Mod UV1 H₂S ≤ 5.000 ppm	2841			H ₂ S	1
xxx = Mod UV1 H₂S ≤ 1%	2856			H ₂ S	1
xxx = Mod UV2	2831			SO ₂ NO ₂ O ₃ Cl ₂ ≤5%	2
xxx = Mod UVRAS	2917			SO ₂ NO ₂ NO	3

Part 4: Casings for Sensor Modules assembled according to Customer Specification			
RMS-Cas-2 ⁽¹⁾	2817	Casing type 2	Suitable for multiple modules 444 x 305 x 145 mm
RMS-Cas-3 ⁽¹⁾	2818	Casing type 3	Suitable for multiple modules 444 x 305 x 189mm
Part 5 Built-in Options			
RMS-O2 ^{(2) (3)}	2795	Oxygen Sensor	0 - 25%
	2767		0 - 100%
RMS-O2-resist ⁽³⁾	2824	Oxygen Sensor H ₂ S resistant	0.5 - 35%
RMS-P ^{(2) (3)}	2771	Pressure Sensor	800 - 1200 mbar abs. Resolution <1 mbar
RMS-P-resist ⁽³⁾	2825	Pressure Sensor H ₂ S resistant	0.2 - 3.5 bar abs. Resolution 2 mbar
RMS-H ^{(2) (3)}	2773	Humidity Sensor	0 - 100% RH
RMS-A/O	2648	Analog Voltage Output 0-2V / 0-5V / 0-10V	4 Analog output ports for 4 separate gas concentrations, 16 bit
RMS-CasHeat	2954	Thermostatted Casing	Heating and thermo- statting of the sensor casing at 50°C

⁽¹⁾ Casing type depends on the type and number of built-in sensor modules

⁽²⁾ Not suitable for SO₂, Cl₂, H₂S

⁽³⁾ Available as supplement to IR or UV sensor only



Casing Type 1



Casing Type 2

Part 6: Accessories / Options

xxx = Cal-ZP-N₂	2805	Calibration Gas N ₂	For zero-point calibration of all gases
xxx = Flow-V	2806	Mini Flow Valve	Control of flow rate for calibration gas bottle incl. manometer
xxx = Cal-CG-Cat1	2948	Calibration with special carrier gas (Ar, H ₂ , He) for gases category 1	Category 1 gases: CO ₂ , CO < 10Vol-%, N ₂ O, CH ₄ , C _n H _m , CF ₄ , SF ₆ , O ₃ , Cl ₂ , NO, NO ₂ , SO ₂ < 10Vol-%
xxx = Cal-CG-Cat2	2949	Calibration with special carrier gas (Ar, H ₂ , He) for gases category 2	Category 2 gases: CO > 10Vol-%, H ₂ S, SO ₂ > 10Vol-%
xxx = Cal-ReCal-Cat1	2950	Recalibration for gases category 1	Category 1 gases: CO ₂ , CO < 10Vol-%, N ₂ O, CH ₄ , C _n H _m , CF ₄ , SF ₆ , O ₃ , Cl ₂ , NO, NO ₂ , SO ₂ < 10Vol-%
xxx = Cal-ReCal-Cat2	2951	Recalibration for gases category 2	Category 2 gases: CO > 10Vol-%, H ₂ S, SO ₂ > 10Vol-%

1.8. Calibrations

Depending on the number of gases to be analysed, a suitable calibration gas must be selected. Furthermore, the type of gas also influences the optimal calibration gas. Because of the almost infinite number of possible combinations, the following list should therefore only be understood as a first overview.

1. One gas: Calibration in N₂

2. Two gases:

a) Binary gas mixture:

Definition: The sum of the concentrations of both gases is 100%.

The Calibration takes place in the respective other gas "gas in gas".
(e.g. CO₂ in CH₄, CH₄ in CO₂)

b) For specific measuring ranges, e.g. CO₂ 0-10 Vol.% and CO 0-20 Vol.%, a calibration of one gas at a time is carried out in nitrogen N₂.

For specific measuring ranges such as CO₂ 0-10 vol.% and CO 0-20 vol.%, a calibration of one gas each is carried out in nitrogen N₂. The calibration would be carried out in the above example as follows:

- CO₂ 10 Vol.% + 90 Vol.% N₂
- CO 20 Vol.% + 80 Vol.% N₂


3. Three and more gases: The cross-sensitivity and, if necessary, the carrier gas dependency of the individual gas components must be checked individually. Then, for example, gas 1 and gas 2 can be defined as binary (and calibrated as "gas in gas") and gas 3 can be calibrated in N₂.

1.9. Measurement of Humid Gases

When measuring humid gases, it is essential to prevent condensation of the water vapour inside the sensor. The following diagram shows the maximum water vapour concentration (= absolute humidity) in volume % as a function of temperature at the standard pressure 1013 mbar.

It can be seen from the diagram:

- At room temperature, no condensation occurs up to a maximum humidity of 2 vol.%.
- Condensation occurs at concentrations above the curve shown.



Condensation can be prevented by the following measures:

- Condensate trap or gas cooler in flow direction in front of the sensor. Please note: This will increase the **dead space** in front of the sensor.
- Installation of a heating element with thermostating in the sensor casing - see para. 7.8. With this heating element, a temperature increase of up to 45°C is possible. Thus, condensation can be prevented up to a water vapour concentration (abs. humidity) of approx. 12 vol.%. **Please note: In addition to heating the sensor, the external gas lines/tubes from the gas source to the sensor unit must also be heated.**

With the optional humidity sensor (see para. b)) the relative humidity can be measured in the measuring range of 0 ~ 100% rH. Furthermore, the absolute humidity can be measured as well in g/m³. Both values are displayed in the supplied software.

2. Disadvantages of electrochemical (EC) sensors in comparison to the optical gas sensors »RITTER MultiGas«

- EC sensors become "blind" over time and then display a constant value, usually zero. This suggests misleadingly a stable zero point.
- EC sensors must therefore be replaced preventatively every 0.5 - 2 years, after replacement the sensor must be recalibrated, as the tolerances of EC sensors are relatively high. This causes supplementary additional costs.
- With EC sensors, mutual influence and deterioration occur by different gases, e.g. NO₂ damages the SO₂ sensor and vice versa.
- EC sensors react very strongly to hydrogen. Such sensors are therefore unusable for accurate measurements even at the smallest concentrations of H₂ in the gas mixture.
- In many countries (e.g. in China) EC sensors are forbidden by law with governmental inspection and approval measurements, because they show too low values if they are contaminated or aged. The user then receives "false positive" values.
- The lifetime of the EC cells is already reduced during storage; therefore, the storage should only be a few weeks.
- The response time (t₉₀) is relatively long compared to the optical measuring methods - mostly about 30 sec. Optical systems are in the range < 5 sec.
- Due to the measuring principle of the EC sensors there is always a chemical reaction between the test gas and the sensor. By this reaction, small quantities of the test gas components are converted. For example, CO is converted into CO₂. With low quantities of test gas, measurements beyond the gas sensor can therefore be influenced because fewer CO molecules are present in the gas sample.

3. Scope of Supply

No. of Items	Item
1	Document folder including ... • Calibration Certificate • Data Sheets • Software Manual
1	Sensor in table top casing
1	Plug-in power supply unit Input: 100-240 V-AC, 50/60 Hz, 0.4 A Output: 24 V-DC, 0.5 A, 12 W or 36 W (H ₂ S / UVRAS sensor)
1	Data acquisition software (on USB memory card)
1	Connection cable sensor ⇒ computer Standard: USB V2.0 A/B, 1 m Option: RS232, D-Sub socket 9pol. on the unit, connection cable [1x D-Sub male 9pin - 1x D-Sub female 9pin], 2 m
3 m	Viton tubing Ø _i 4 mm / Ø _o 6 mm

4. Setup

1. Unpack all items carefully
2. Install the data acquisition software according to the "RITTER MultiGas Software Manual".
Don't start the software at this point in time.
3. Positioning of the sensor:

- a) **The sensor must not be used in explosive environments** or in harsh environmental conditions (e.g. high condensing humidity, strong air currents, aggressive atmospheres, outdoors without housing).
- b) **If the sensor is used in combination of a RITTER MilliGascounter, the sensor must be positioned (in flow direction) in front of the MilliGascounter if this MilliGascounter is filled with HCl solution. Otherwise, the sensor could be damaged by the HCl vapour coming out of the MilliGascounter.**

- c) Place the sensor next to the gas source
Please note: A tube connection between gas source and sensors that is as short as possible allows a small dead space created by the tubing. In turn, a small dead space enables a fast response time of the sensor.
4. Connect the power supply unit to the mains and to the »DC 24 V« socket on the rear (casing type 1) or front (casing type 2) of the sensor.
5. Connect the data acquisition cable to the respective socket at the rear side of the sensor and to the computer:
 - a) USB cable to the socket "USB"
 - b) RS232 cable to the socket "RS 232" (option)
6. Connect the gas source to the gas inlet port of the sensor by using the provided tubing as follows:
 - a) Remove the screw caps as well as the red rubber caps of the transport lock from the gas inlet and gas outlet.
 - b) Slide the screw cap onto the end of the gas inlet tube with the thread of the screw cap facing the tube end.
 - c) Push the tube onto the cone in the center of the gas inlet port.
 - d) Slide the screw cap to the gas inlet port and screw it hand-tightly in place.
7. If applicable: Connect the gas outlet port of the sensor to other components such as gas sampling bags, exhaust tubing etc.
 - The connection of the tube to the gas outlet port takes place in the same way as described above.
8. Switch-on the power switch at the front side of the sensor casing.
9. Start the software and open the COM port(s) of the connected sensor module(s) according to the software manual.

The »RITTER MultiGas« Sensor is now ready for use.

5. Infrared Sensor


5.1. Description

Gas analysis based on the NDIR technique is an established method to determine the concentrations of gases in complex mixtures. The »RITTER MultiGas« sensors use new optical components for optimal analysis results: Up to 3 optical filters analyse the gas which flows through the sensor as one gas stream. The optional oxygen, pressure, and humidity sensors are in-line with the same single gas stream.

The individual internal modules are sealed by means of O-ring connections.

In order to achieve an optimum adaptation to the required measuring range, the lengths of the modular measurement cells (= cuvettes) can be implemented in the range of 5 mm (large measurement range in percentage level) up to 250 mm (small measurement range in ppm level).

Cuvettes with a length ≥ 20 mm are coated with a resistant gold layer in order to improve the reflection properties for



low concentration level detection. Cuvettes used with aggressive gases are gold-coated as well.

The other internal mechanical parts are made out of aluminium, optionally out of stainless steel.

For fast response applications the measuring system delivers a stable result within $t_{90} \approx 3$ seconds.

The entire unit can be disassembled for easy maintenance/service.

5.2. Applications

- Biogas / natural gas analysis
- Environmental and Process Measurement
- TOC analysers
- Continuous Emission Monitoring (CEM)
- Elemental analysis
- Industrial gas analysis

5.3. Specifications

General Features	
Measurement technology	Innovative NDIR Sensor (non-dispersive infrared sensor)
Detectable gases	CO ₂ , CO, N ₂ O, NO, CH ₄ , C _n H _m *, CF ₄ , SF ₆
Number of simultaneously detectable gases	max. 3 per sensor unit
Measurement ranges	See par. 5.4 or https://www.ritter.de/en/products/sensors#ranges
Flow rate range	5 ~ 300 ltr/h For higher flow rates the sensor can be operated in bypass
Max. gas inlet pressure	300 mbar
Pressure loss (without additional optional sensors)	10 @ 100 / 35 @ 200 / 70 @ 300 [mbar @ ltr/h]
Temperature compensation	Yes

Data acquisition software	Yes
Lifetime of IR radiation source	> 40 000 h
Measurement cuvette	Aluminium, with measurement ranges $\leq 1\%$ gold-plated inside
Cuvette sealing	Viton O-ring
Internal tubing	FKM / Viton TM (fluorinated rubber)
Casing	High-quality table-top casing, aluminium
Dimensions	W x H x L 171 x 86 x 290 mm
Weight	approx. 2 kg
Gas connections	PVDF screw-type tube connection for tube $\varnothing_i$ 4mm, $\varnothing_o$ 6 mm

Measurement Data

Linearity error	< $\pm 1\%$ F.S.
Repeatability	$\pm 0.5\%$ F.S.
Long term stability zero	< $\pm 2\%$ F.S. / week
Long term stability span	< $\pm 2\%$ F.S. / month
Temperature influence of zero point	< 1% F.S. / 10K
Temperature influence of span	< 2% F.S. / 10K
Cross sensitivity	< 2% F.S.
Pressure influence	< 1.5% / 10hPa of reading
Warm-up time	2 min
Response time (t_{90})	≈ 3 sec
Sampling frequency by software	≤ 10 Hz
Detection limit ($3 \cdot \sigma$)	See par. 5.5 or https://www.ritter.de/en/products/sensors#limits
Resolution	0.5 x detection limit
Water vapour	No influence with measurements of CO ₂ and CH ₄

Electrical Features

Power supply	24 VDC, incl. power plug 100~240 VAC, 50/60 Hz $\Rightarrow$ 24 VDC
Average power consumption	< 1W
Interfaces	standard: USB options: RS232, CANbus, CANopen incl. data transmission cable 1 m
Analogue voltage output (option)	0-2V / 0-5V / 0-10V

Climatic Conditions

Operating temperature +15 ~ +45 °C

Storage temperature -20 ~ +60 °C

Operating pressure 800 ~ 1200 hPa (mbar)

Ambient humidity 0 ~ 95% rel. humidity
Condensing inside of sensor must be prevented!

* Analysis of C_nH_m :

The calibration of sensors for C_nH_m will be performed with propane. Aromatic hydrocarbons are also measured but with a different weighting. This means that the sensitivity of the sensor is significantly smaller with these gases than with other hydrocarbons.

5.4. Table of Standard Measurement Ranges

Full Scale ¹⁾	CO ₂	CO	N ₂ O	CH ₄	C _n H _m ²⁾	CF ₄	SF ₆	H ₂ O
100 vol. %	x	x	x	x	x	x	x	
50 vol. %	x	x	x	x	x		x	
30 vol. %		x	x	x	x		x	
20 vol. %	x							
10 vol. %	x	x		x	x			
5 vol. %	x	x		x	x			
2 vol. %								x
1 vol. %	x	x		x	x			
5,000 ppm	x	x		x	x		x	
2,000 ppm	x	x	x	x	x			
1,000 ppm	x	x	x	x	x		x	
500 ppm	x	x	x					
300 ppm			x					
100 ppm	x		x				x	
50 ppm	x						x	

¹⁾ Measurement range = Zero to Full Scale (FS)

Sensors are calibrated in the same range "Zero to Full Scale (FS)"

²⁾ Calibration with Propane

Other ranges on request

5.5. Table of Detection Limits (= 3 σ) in Percent of Full Scale (FS)

Full Scale	CO ₂	CO	N ₂ O	CH ₄	C _n H _m	CF ₄	SF ₆	H ₂ O
100 vol. %	< 0.1%	< 0.2%	< 0.1%	< 0.1%	< 0.1%	< 0.2%	< 0.1%	
50 vol. %	< 0.1%	< 0.2%	< 0.1%	< 0.1%	< 0.1%	< 0.2%	< 0.1%	
30 vol. %		< 0.2%	< 0.1%	< 0.1%	< 0.1%	< 0.2%	< 0.1%	
20 vol. %	< 0.1%							
10 vol. %	< 0.1%	< 0.2%		< 0.1%	< 0.2%			
5 vol. %	< 0.1%	< 0.2%		< 0.1%	< 0.2%			
2 vol. %								0,005%
1 vol. %	< 0.1%	< 0.2%		< 0.1%	< 0.2%			
5,000 ppm	< 0.1%	< 0.2%		< 0.1%	< 0.2%			
2,000 ppm	< 0.1%	< 0.3%	< 0.1%	< 0.3%	< 0.5%			
1,000 ppm	< 0.1%	< 0.5%	< 0.1%	< 0.5%	< 0.5%			
500 ppm	< 0.1%	< 0.5%	< 0.1%					
300 ppm	< 0.1%		< 0.1%					
100 ppm	< 0.3%		< 0.3%					
50 ppm	< 0.3%							

Definition of Detection Limit

The Detection Limit is the smallest measurement value which can be obtained with a specific uncertainty. This uncertainty includes the resolution, noise and stability of the gas sensor for a specific gas and specific measurement range. For evaluation of the detection limit value, several single measurements are taken at the identical measurement conditions. With the obtained single measurement results the standard deviation "Sigma" (σ) is calculated. The values given in the table equal the triple amount of Sigma.

5.6. Recalibrations

The following recalibration intervals are recommended for IR sensors:

- Zero-point: Weekly with air or an inert gas, e.g. Nitrogen
The recalibration of the zero point is described in the software manual.
- End-point (full scale): Every 3 months with suitable calibration gas

6. Ultraviolet Sensors

6.1. Description


The »RITTER MultiGas« UV sensor is the world's first gas analysis module based on miniaturized UV-LEDs. The stability and lifetime of these UV-LEDs enables high-precision gas analyses down to the ppm range. By using two UV-LEDs two gases can be detected simultaneously. Furthermore, with this approach measuring ranges from ppm to Vol.-% can be realized.

In the spectral range from 200 nm to 500 nm, nitrogen oxides (NO+NO₂), aromatic hydrocarbons, hydrogen sulphide, ozone, sulphur dioxide and chlorine can be reliably detected with this new sensor platform.

The entire unit can be disassembled for easy maintenance/service.

The individual internal modules are sealed by means of O-ring connections.

In order to achieve an optimum adaptation to the required measuring range, the lengths of the modular measurement cells (= cuvettes) can be implemented in the range of 5 mm (large measurement range



in percentage level) up to 250 mm (small measurement range in ppm level).

Cuvettes with a length ≥ 20 mm are coated with a resistant gold layer in order to improve the reflection properties for low concentration level detection. Cuvettes used with aggressive gases are gold-coated as well.

The internal mechanical parts are made out of aluminium, optionally out of stainless steel.

For fast response applications the measuring system delivers a stable result within $t_{90} \approx 1-2$ seconds.

6.2. Applications

- Biogas / natural gas analysis
- Environmental and Process Measurement
- TOC analysers
- Continuous Emission Monitoring (CEM)
- Elemental analysis
- Industrial gas analysis

6.3. Specifications NDUV Sensor

General Features

Measurement technology	Innovative NDUV Sensor (non-dispersive ultraviolet sensor)
Detectable gases	SO ₂ , NO ₂ , O ₃ , Cl ₂
Number of simultaneously detectable gases	max. 2
Measurement ranges	See par. 6.6 or https://www.ritter.de/en/products/sensors#ranges
Flow rate range	5 ~ 300 ltr/h For higher flow rates the sensor can be operated in bypass
Max. gas inlet pressure	300 mbar

Pressure loss (without additional optional sensors)	10 @ 100 / 35 @ 200 / 70 @ 300 [mbar @ ltr/h]
Temperature compensation	Yes
Data acquisition software	Yes
Lifetime of UV radiation source	> 8 000 h
Measurement cuvette	Stainless steel with silicone coating inside
Cuvette sealing	Viton O-ring
Internal tubing	FKM / Viton TM (fluorinated rubber)
Casing	High-quality table-top casing, aluminium
Dimensions	W x H x L 464 x 189 x 305 mm
Weight	approx. 6.5 ⁺ kg
Gas connections	PVDF screw-type tube connection for tube Ø _i 4mm, Ø _e 6 mm

Measurement Data

Linearity error	< ± 1% F.S.
Repeatability	± 0.5% F.S.
Long term stability zero N ₂	< ± 1% F.S. / 24h
Long term stability span	< ± 1% F.S. / month
Temperature influence of zero point	< 1% F.S. / 10K
Temperature influence of span	< 2% F.S. / 10K
Cross sensitivity	< 2% F.S.
Pressure influence	< 1.5% / 10hPa of reading
Warm-up time	1 min (initial), <60 min for full specification
Response time (t ₉₀)	1.5 - 15 sec
Sampling frequency by software	≤ 10 Hz
Detection limit (3·σ)	See par. 6.7 or https://www.ritter.de/en/products/sensors#limits
Resolution	0.5 x detection limit

Electrical Features

Power supply	24 VDC, incl. power plug 100~240 VAC, 50/60 Hz ⇒ 24 VDC
Supply current (peak)	<0.4 A
Average power consumption	< 7.5 W
Interfaces	standard: USB options: RS232, CANbus, CANopen incl. data transmission cable 1 m
Analogue voltage output (option)	0-2V / 0-5V / 0-10V


Climatic Conditions

Operating temperature	+25 ~ +45 °C
Storage temperature	–20 ~ +60 °C
Operating pressure	800 ~ 1200 hPa (mbar)
Ambient humidity	0 ~ 95% rel. humidity Condensing inside of sensor must be prevented!

6.4. Specifications H₂S Sensor

The H₂S sensor works according to the principle of non-dispersive UV absorption (NDUV).

The measuring wavelength used is in the lower nanometer range, which enables interference-free measurements with water vapour and hydrocarbons. This makes the sensor ideal for use in biogas measurements, as accompanying gases such as NH₃, CO₂, CH₄ and H₂O do not affect the accuracy of the measurements.



General Features

Measurement technology	Innovative NDUV Sensor (non-dispersive ultraviolet sensor) Dual beam
Detectable gases	H ₂ S, SO ₂ , CH ₄
Number of simultaneously detectable gases	max. 2
Measurement ranges	See par. 6.6 or https://www.ritter.de/en/products/sensors#ranges 5 ~ 300 ltr/h
Flow rate range	For higher flow rates the sensor can be operated in bypass
Max. gas inlet pressure	300 mbar
Pressure loss (without additional optional sensors)	10 @ 100 / 35 @ 200 / 70 @ 300 [mbar @ ltr/h]
Temperature compensation	Yes
Data acquisition software	Yes
Lifetime of UV radiation source	> 8 000 h
Measurement cuvette	Stainless steel with silicone coating inside
Cuvette sealing	Viton O-ring
Internal tubing	FKM / Viton™ (fluorinated rubber)
Casing	High-quality table-top casing, aluminium
Dimensions of module	W x H x L 300 x 100 x 81 mm
Weight of module	approx. 1670 g

Gas connections	PVDF screw-type tube connection for tube $\varnothing_i$ 4mm, $\varnothing_o$ 6 mm
-----------------	---

Measurement Data

Linearity error	$< \pm 1\%$ F.S.
Repeatability	$\pm 0.5\%$ F.S.
Long term stability zero N ₂	$< \pm 1\%$ F.S. / 24h
Long term stability span	$< \pm 1\%$ F.S. / month
Temperature influence of zero point	$< 1\%$ F.S. / 10K
Temperature influence of span	$< 2\%$ F.S. / 10K
Cross sensitivity	$< 2\%$ F.S.
Pressure influence	$< 1.5\%$ / 10 hPa of reading
Warm-up time	1 min (initial), < 15 min for full specification
Response time (t ₉₀)	1.5 - 15 sec
Sampling frequency by software	≤ 10 Hz
Detection limit (3· σ)	See par. 6.7 or https://www.ritter.de/en/products/sensors#limits
Resolution	0.5 x detection limit

Electrical Features

Power supply	24 VDC, incl. power plug 100 ~ 240 VAC, 50/60Hz $\Rightarrow$ 24 VDC
Supply current (peak)	< 0.4 A
Average power consumption	< 7.5 W
Interfaces	standard: USB options: RS232, CANbus, CANopen incl. data transmission cable 1 m
Analogue voltage output (option)	0-2V / 0-5V / 0-10V

Climatic Conditions

Operating temperature	+5 ~ +45 °C
Storage temperature	-20 ~ +60 °C
Operating pressure	800 ~ 1200 hPa (mbar)
Ambient humidity	0 ~ 95% rel. humidity Condensing inside of sensor must be prevented!

6.5. Specifications UVRAS Sensor

For the detection of NO an EDL (electrodeless gas discharge lamp) is used. In the EDL, N₂ and O₂ are converted to NO and produce a selective UV radiation. With this radiation, a cross-sensitivity-free NO measurement is made possible. This method is called UV resonance absorption spectroscopy (UVRAS).

A combination of both the UVRAS and NDUV technology allows the simultaneous gas analysis of NO, NO₂ and SO₂ in the lower ppm range which is particularly important in flue gas analysis (Continuous Emission Monitoring, CEM).



General Features

Measurement technology	UV resonance absorption spectroscopy (UVRAS)
Detectable gases	NO, NO ₂ , SO ₂
Number of simultaneously detectable gases	max. 3
Measurement ranges	See par. 6.6 or https://www.ritter.de/en/products/sensors#ranges
Flow rate range	5 ~ 300 ltr/h For higher flow rates the sensor can be operated in bypass
Pressure loss (without additional optional sensors)	10 @ 100 / 35 @ 200 / 70 @ 300 [mbar @ ltr/h]
Additional heated and thermostated enclosure for UVRAS module inside of sensor casing	Yes
Temperature compensation	Yes
Data acquisition software	Yes
Lifetime of UV radiation sources	LED > 20 000 h (NO ₂ , SO ₂) EDL > 8 000 h (NO)
Measurement cuvette	Stainless steel with silicone coating inside
Cuvette sealing	Viton O-ring
Internal tubing	FKM / Viton™ (fluorinated rubber)
Casing	High-quality table-top casing type 2, aluminium
Dimensions	W x H x L 464 x 189 x 305 mm
Weight	approx. 6.5 ⁺ kg
Gas connections	PVDF screw-type tube connection for tube Ø _i 4mm, Ø _o 6 mm

Measurement Data

Linearity error	< ± 1% F.S.
Repeatability	± 0.5% F.S.
Long term stability zero N ₂	< 3 ppm / 24h
Long term stability span	< ± 1% F.S. / month
Temperature influence of zero point	< 1% F.S. / 10K
Temperature influence of span	< 2% F.S. / 10K
Cross sensitivities	500 ppm NO ₂ < 2ppm 100 ppm SO ₂ < 2ppm 100 ppm N ₂ O < 10ppm 20°C D.P. H ₂ O < 10ppm
Pressure influence	< 1.5% / 10hPa of reading
Warm-up time	1 min (initial), <60 min for full specification
Response time (t ₉₀)	1.5 ~ 15 sec
Sampling frequency by software	≤ 10 Hz
Detection limit (3·σ)	See par. 6.7 or https://www.ritter.de/en/products/sensors#limits
Resolution	0.5 x detection limit

Electrical Features

Power supply	24 VDC, incl. power plug 100~240 VAC, 50/60 Hz ⇒ 24 VDC
Supply current (peak)	1.5 A
Inrush current	0.2 ~ 0.7 A
Power consumption (peak)	36 W
Interfaces	standard: USB options: RS232, CANbus, CANopen incl. data transmission cable 1 m
Analogue voltage output (option)	0-2V / 0-5V / 0-10V

Climatic Conditions

Operating temperature	+5 ~ +40 °C
Storage temperature	-20 ~ +60 °C
Operating pressure	800 ~ 1200 hPa (mbar)
Ambient humidity	0 ~ 95% rel. humidity Condensing inside of sensor must be prevented!

6.6. Table of Standard Measurement Ranges

Full Scale ¹⁾	O ₃	CL ₂	H ₂ S	NO	NO ₂	SO ₂
30 vol. %		x				
10 vol. %		x				x
5 vol. %		x				x
1 vol. %		x	x			x
5,000 ppm		x	x	x	x	x
2,000 ppm	x	x	x	x	x	x
1,000 ppm	x	x	x	x	x	x
500 ppm	x	x	x	x	x	x
300 ppm				x	x	x
100 ppm	x	x	x		x	x
50 ppm	x				x	x
10 ppm	x				x	x
1 ppm	x					

¹⁾ Measurement range = Zero to Full Scale (FS)
Sensors are calibrated in the same range "Zero to Full Scale (FS)"
Other ranges on request

6.7. Table of Detection Limits (= 3 σ) in Percent of Full Scale

Full Scale	O ₃	CL ₂	H ₂ S	NO	NO ₂	SO ₂
100 vol. %						
50 vol. %						
30 vol. %		< 0.1%				
20 vol. %						
10 vol. %		< 0.1%				< 0.1%
5 vol. %		< 0.1%				< 0.1%
1 vol. %			< 0.1%			
5,000 ppm			< 0.1%			
2,000 ppm	< 0.1%	< 0.1%	< 0.1%	< 0.1%	< 0.1%	< 0.1%
1,000 ppm	< 0.1%	< 0.1%	< 0.2%	< 0.1%	< 0.1%	< 0.1%
500 ppm	< 0.2%	< 0.2%	< 0.3%	< 0.2%	< 0.2%	< 0.1%
300 ppm				< 0.2%	< 0.2%	< 0.1%
100 ppm	< 0.5%		< 0.5%		< 0.5%	< 0.3%
50 ppm	< 0.5%				< 0.5%	< 0.3%
10 ppm	< 0.5%				< 0.5%	< 0.3%

Definition of Detection Limit

The Detection Limit is the smallest measurement value which can be obtained with a specific uncertainty. This uncertainty includes the resolution, noise and stability of the gas sensor for a specific gas and specific measurement range. For evaluation of the detection limit value,

several single measurements are taken at the identical measurement conditions. With the obtained single measurement results the standard deviation "Sigma" (σ) is calculated. The values given in the table equal the triple amount of Sigma.

6.8. Recalibrations

The following recalibration intervals are recommended for UV sensors:

- Zero-point:
 - Concentrations < 300 ppm: Every 48 hours with inert gas, e.g. Nitrogen
 - Concentrations $\geq$ 300 ppm: Every 24 hours with inert gas, e.g. NitrogenThe recalibration of the zero point is described in the software manual.
- End-point (full scale): Every 3 months with suitable calibration gas

7. Options (installed inside of the sensor casing)

Please note: All optional sensors are available in addition to a RITTER »MultiGas« NDIR or NDUV sensor only.

7.1. Oxygen Sensor

The measured oxygen concentration is displayed in the provided software. The sensor is built into the casing of the RITTER MultiGas sensors.

The following versions are available:

- a) For non-aggressive gases
- b) For H₂S and similar acid gases



Specifications

Dimensions [H x W x L] 65.4 mm x 31.7 mm x 56.6 mm

Weight 70 g

Tube connector 4/6 mm tube

Measurement range	0 – 25 Vol.% O ₂	0 – 100 Vol.% O ₂	0.5 – 35 Vol.% O ₂
Application	Biogas, Automotive exhaust gas analyser	Industrial, fully CO ₂ resistant	Industrial, fully CO ₂ resistant, shows high resistance to acid gases
Medium contact materials	ABS, FKM, PPS, PTFE, stainless steel	ABS, PVC, PPS, PTFE, stainless steel	ABS, PVC, PPS, PTFE, stainless steel
Expected operating life	> 1,000,000 Vol.% O ₂ h	~ 1,200,000 Vol.% O ₂ h	~ 1,200,000 Vol.% O ₂ h
Sensor lifetime	4 years @ ambient air, depending on application	6 years @ ambient air, depending on application	6 years @ ambient air

Measurement Data¹

Resolution	0.1 Vol.%	0.1 Vol.%	0.1 Vol.%
Response time (t ₉₀)	< 3.5 s	< 10 s	< 5 s
Drift ²	< 1% per month	< 1% per month	< 3% per month

¹ related to Pa = 1013 hPa; Ta = 25°C; RH = 50%; flow = 2.5l/min

² averaged across 12 months

Linearity Error	0 - 2 Vol.% O ₂ : ± 0.1 abs. 2.1 - 100 Vol.% O ₂ : ± 0.05 rel.	0 - 2 Vol.% O ₂ : ± 0.1 abs. 2.1 - 35 Vol.% O ₂ : ± 0.05 rel.
-----------------	---	--

Repeatability ³	± 1 Vol.% O ₂	± 1 Vol.% O ₂
----------------------------	--------------------------	--------------------------

Influence of Humidity	-0.03 % rel. O ₂ reading / % RH	-0.03 % rel. O ₂ reading / % RH	-0.03 % rel. O ₂ reading / % RH
-----------------------	--	--	--

Interferences	CO ₂ : up to 20 Vol.% CO: up to 2000 ppm NO _x : up to 5000 ppm HC: up to 5000 ppm N ₂ O: up to 500 ppm	< 20 ppm O ₂ response to: 100 Vol.% CO 100 Vol.% CO ₂ 100 Vol.% C ₃ H ₈ 3000 ppm NO in N ₂ 1000 ppm C ₆ H ₆ in N ₂ 500 ppm SO ₂ in N ₂ < 100 ppm O ₂ response to: 3000 ppm C ₂ H ₆ O 3000 ppm C ₄ H ₁₀ S < 200 ppm O ₂ response to: 3000 ppm C ₂ H ₆ S ₂ < 400 ppm O ₂ response to: 100 Vol.% H ₂ < 500 ppm O ₂ response to: 2000 ppm H ₂ S in N ₂	< 20 ppm O ₂ response to: 100 Vol.% CO 100 Vol.% CO ₂ 100 Vol.% C ₃ H ₈ 1000 ppm C ₆ H ₆ in N ₂ 2000 ppm H ₂ S in N ₂ < 20000 ppm O ₂ response to: 3000 ppm NO in N ₂ 1000 ppm H ₂ in N ₂ 500 ppm SO ₂ in N ₂
---------------	---	---	---

Climatic Conditions

Operating temperature	0 – 40 °C intermittent 40 – 50 °C -20 – 40 °C	0 – 45 °C	0 – 50 °C
-----------------------	---	-----------	-----------

Storage temperature	5 – 25 °C recommended 40 – 50 °C max. 1 week	5 – 30 °C recommended -20 – 50 °C maximum
---------------------	---	--

Air pressure	650 – 1250 hPa (mbar)	700 – 1250hPa (mbar)	600 – 1250 hPa (mbar)
--------------	-----------------------	----------------------	-----------------------

Ambient humidity	0 – 95% rel. humidity (not condensing)
------------------	--

³ @ 100 Vol.% O₂ applied for 5 min

7.2. Pressure Sensor

A change in the gas and/or atmospheric pressure causes a change in the number of molecules per volume and thus a change in gas density. This density change in turn has a significant influence on the result of the concentration measurement by the sensor. By measuring the gas pressure inside of the sample cell (cuvette), the value of the concentration measurement is compensated / corrected.

The pressure sensor enables an **increase in the accuracy** of the gas analysis measurement by a factor of 15:

- Without pressure sensor: **±1.5%** per pressure difference of 10 hPa
- With pressure sensor: **±0.1%** per pressure difference of 10 hPa


The pressure sensor is built into the casing of the RITTER MultiGas sensors. The measured pressure value is displayed in the provided software.

Pressure Sensor for IR Sensors

The pressure sensor for IR sensors is provided as a **standard equipment with all IR sensors**.

Specifications

- Measuring range: 800 ~ 1,200 mbar abs.
- Measurement accuracy ±1% of span (full scale)
- Resolution: <1 mbar
- Response time (t_{90}): 1 s
- Incl. temperature compensation
- For non-aggressive gases



Pressure Sensor for UV Sensors

The pressure sensor for UV sensors is provided as an **option on request**.

- a) For non-aggressive gases

Specifications

- Measuring range: 600 ~ 1.100 mbar abs.
- Measurement accuracy ±1% of span (full scale)
- Resolution: 0.1 mbar
- Response time (t_{90}): 1.5 ~ 15 s
- Incl. temperature compensation

- b) For H₂S and similar acid gases


Specifications

- Measuring range: 0.2 ~ 3.5 bar abs.
- Measurement accuracy ±1% of span (full scale)
- Resolution: 2 mbar
- Response time (t_{90}): 1.5 ~ 15 s
- Incl. temperature compensation



7.3. Humidity Sensor

The humidity sensor is a sensor module available as option in addition to a RITTER MultiGas NDIR or NDUV sensor. The sensor is built into the casing of the RITTER MultiGas sensors. The measured humidity values (absolute and relative) are displayed in the provided software.



Specifications


- Polymer humidity sensor
- Measuring range: 0-100% rH
- Measurement accuracy $\pm 2\%$ rH of span (full scale)
- Resolution: $\pm 1\%$ RH
- Response time (t_{90}): 12 s
- Incl. temperature compensation
- Indicated values (in software): absolute [% absH] and relative humidity [% rH]

Display of the measured values of the humidity sensor in the data acquisition software "MARS" and their output as analog voltage signal:

The following measured values of the humidity sensor are displayed in the software in the data logging window (*Data Logging Tab*) as follows:

- Absolute humidity [g/m^3]
- Relative humidity [%]
- Concentration [ppm]
- Temperature in humidity sensor [$^{\circ}\text{C}$]

From these values, the measured **concentration** is output as an analog voltage signal on channel 4 (grey wire) for NDIR sensors, and on channel 2 (green wire) for NDUV sensors. (See also table above)



The output voltage value is determined as follows:

In the software, the measuring range for the concentration value is preset to 30,000 ppm in the Channel 4 configuration window (*Config Tab CH4*). (Can be changed individually)

The measured value of the concentration in the data acquisition window is set in relation to this preset value. This quotient is multiplied by the voltage level of the analog output and output as a voltage value.

Example:


Measured value: 12,000 ppm

Measuring range in *Config Tab CH4*:
30,000 ppm

Quotient $12,000 / 30,000 = 0.4$

Voltage level of the analog output: 0 - 10V

Output voltage: $0.4 \times 10V = 4.0V$



MARS - Multi Analyse Research Software

Interface
☒ RS232 ☐ COM Port ☐ CAN ☐ COM10 ☐ CAN ☐ COM10

Serial Number: Detector connected: ?

Data Logging Basic configuration Config CH1 Config CH2 Config CH3 **Config CH4** Production S

Linearisation_Ch4 (x*) 0.00000 0.00000 0.00000 0.00000
 0.00000 0.00000

Enable Ln-Function Ch4 ☐

TempCompNullpoint_Ch4 (b) 0.00000 0.00000 0.00000 10E-6 ☐

TempCompEndpoint_Ch4 (b) 0.00000 0.00000 0.00000 10E-6 ☐

Merge CH3 to CH4 ☐

Endpoint Calibration Ch4

Endpoint Factor Ch4 1.00000


µg Factor

Measurement Unit Ch4 ppm

Measurement Range Ch4 30000.00000

Zero Detector Ch4

7.4. Automatic Zero-point Adjustment by Air




Function:

- Automatic time-controlled zero-point adjustment for long-term stable measurements
- Setting of time interval in supplied data acquisition software
- Zero setting with ambient air
- Including pressure compensation


Specifications:

- Time interval for zero-point adjustment: from 1 hour to 1 month
- Response time of zero setting: approx. 15 sec
- Pressure compensation included
- Pressure drop: 25 mbar at 1 ltr/min
- Flow rate of zero gas (ambient air) by inbuilt mini pump: 0.35 ltr/min
- Medium contact materials: EPDM / PPS / Al Mg 4.5 Mn (anodized)
- Power consumption: 2.2 W
- Operating temperature: 5°C ~ 55°C
- Optional module, installed in the casing of the »RITTER MultiGas« sensor
- Operating temperature: 5°C ~ 55°C
- Built-in as supplement in table-top casing of gas sensor »RITTER MultiGas«

Pressure Drop



Circuit Diagram



7.5. Automatic Zero-point Adjustment by Nitrogen




Function:

- Automatic time-controlled zero-point adjustment for long-term stable measurements
- Setting of time interval in supplied data acquisition software
- Zero setting with ambient air
- Including pressure compensation


Specifications:

- Time interval for zero-point adjustment: from 1 hour to 1 month
- Response time of zero setting: approx. 15 sec
- Pressure compensation included
- Pressure drop: 25 mbar at 1 ltr/min
- Flow rate of zero gas (ambient air) by inbuilt mini pump: 0.35 ltr/min
- Medium contact materials: EPDM / PPS / Al Mg 4.5 Mn (anodized)
- Power consumption: 2.2 W
- Operating temperature: 5°C ~ 55°C
- Optional module, installed in the casing of the »RITTER MultiGas« sensor
- Operating temperature: 5°C ~ 55°C
- Built-in as supplement in table-top casing of gas sensor »RITTER MultiGas«

Pressure Drop



Circuit Diagram



7.6. LCD-Display for Sensor Casing

The enables the display of the most important sensor measurement values:

- 4 gas concentrations
- Optional: Oxygen concentration
- Gas temperature
- Gas pressure
- Relative / absolute humidity

7.7. Analog Voltage Output

The Analog Voltage Output Module allows the simultaneous output of up to 4 different measurement values at 4 separate channels:

- 4 parallel signal outputs
- The measured values of the oxygen and humidity sensor can only be transmitted alternatively
- Resolution: 16 bit
- Update rate: 1 sec


For connection to an analog data acquisition device the analog voltage output port provides alternatively the following voltage levels:

- 0 – 2V
- 0 – 5V
- 0 – 10V


The voltage range is pre-set according to the order and cannot be changed by the user.

Please note the assignments of the respective gas concentration measurements to the respective channels. These assignments are stated in the Calibration Protocol of the sensor unit.

Assignment of contacts to the channels of the gas concentration measurements at the



socket and plug of the Analog Voltage Output:
View to socket



View to plug



Contact No.	Assignment to channel	Wire Colors of provided Connection Cable
1	channel 1	white
2	channel 2	green
3	channel 3	yellow
4	channel 4	grey
5	ground	brown

Analogue output signal of **oxygen and humidity sensors**:
(Only one of these two signals can be transmitted alternatively.)


- NDIR sensors: Channel 4 (grey wire)
- NDUV sensors: Channel 2 (green wire)

The connection cable is delivered with open wires to be connected to the analog data acquisition device of the user. The wire colours are listed in the above table.

7.8. Heating and Thermostating of Sensor Casing

In order to avoid condensation of humid gas inside of the sensor, the sensor and sensor casing can be heated and thermostatted. As standard, the temperature controller is pre-set to 45°C; lower temperatures are possible on request.

The heating element and temperature controller are mounted at the support of the measurement cuvette.



Features:

- Temperature controller 30 ~ 45°C
- Control accuracy ± 0.2 K
- Heating power 12 Watt

7.9. Tube Port for Flushing of Casing

Optional additional tube connection at the sensor casing for flushing the casing interior

Functions:

- The gas lines inside the casing and the measuring cuvette are gas-tight by means of O-rings and other gas-tight connections. However, as with all connections, a micro leakage rate cannot be ruled out. With the additional tube connection, an **accumulation of the measuring gas inside the casing can be prevented**.
- In case of **toxic and/or aggressive gases**, a suction line can be connected via this tube connection creating an under-pressure inside the casing and thus preventing the gas from escaping from the casing.
- In case of **ignitable gases** (methane, hydrogen, etc.), a pressure or suction line can be connected via this tube connection creating an over- or under-pressure inside the casing thus preventing the accumulation of an ignitable gas mixture inside the casing.

Notes:

- If a **sensor and/or casing heating is used at the same time**, the interior temperature is lowered, which may result in condensation of the gas to be measured.
- H₂S can be perceived as an unpleasant smell even in the ppb range.

Design:

- 1 x screw-type tube connection for connection of tube Øi 4 mm, Øo 6 mm, open to the inside

7.10. Filter

The filter for gas particles is provided as an accessory for those applications where no filtration of the gas is done on the customer side. It is designed as an inline filter which can easily be inserted in the gas tubing in front of all RITTER »MultiGas« Sensors.



Please note:

- It is essential that particles do not enter the sensor as this will block the gas channels inside the sensor.
- Please ensure that this filter is suitable for your application in terms of material resistance and filtration.

Therefore, one filter per gas line is delivered with each sensor. Spare filters can be ordered from your national RITTER partner – see <https://www.ritter.de/en/worldwide/>

Application

- Particle filter
- Accessory for RITTER »MultiGas« Sensors


Features & Benefits

- Low pressure drop

- Gas connections for 4/6 mm tube
- Filtration of particles > 0.1 µm up to 99.9998%

Specifications

Dimensions	Ø 26.4 mm, overall length 92 mm
Weight	26 g
Tube connector	4/6 mm tube
Flowrate	0 ~ 4 ltr/min
Application	Particle filter
Housing volume	110 ml
Medium contact materials	Polyamide (housing), stainless steel (connector), FKM (O-ring)
Expected operating life	Depends on gas contamination
Pressure drop (mbar) @ 1 ltr/min	<6 mbar
Operating temperature	-5 °C ~ 60 °C
Storage temperature	-20 °C ~ 110 °C
Max. pressure	8 bar



8. Display for Sensor Casing Type 1

Sensor casings of Type 1 can optionally be equipped with an LCD display. The display shows:

- Measurement values from the IR, UV, and optional sensor modules
- Temperature of the IR and UV sensor modules
- Gas pressure



8.1. Stand Bracket / Carrying Handle

The stand bracket, which also serves as a carrying handle for transporting the sensor, allows the casing to be positioned so that the display can be viewed at an optimal angle for the user.

To adjust the desired tilt angle, the locking mechanism of the bracket joints can be released by simultaneously pressing the round release buttons inward – see red arrows in the photo above.

Note: It is possible that after positioning the bracket may not lie exactly parallel to the surface. This is due to unavoidable clearance in the bracket's swivel joint and can be corrected by gently twisting the two bracket arms against each other.


8.2. Startup Phase

After switching on the sensor's main power switch, the display module initializes automatically during the startup phase and initially shows the welcome screen.

Next, the current software version is shown, which defines the respective range of functions available.

In the third phase of initialization, the external components are initialized. At this point, among others the CANopen bus and the sensor modules are connected. If any of the components encounters an unexpected internal error, the initialization process will pause at that point until a connection can be successfully established.

The successful connection of all components is confirmed by this screen:



8.3. Measurement Display Screen

The measurement display is the main screen of the display. Measurement values are shown in four lines at a time. The »Up« and »Down« buttons can be used to scroll through the individual lines.

The following values are displayed on this screen:

- Ch 1: Channel 1
- Ch 2: Channel 2
- Ch 3: Channel 3
- Ch 4: Channel 4
- Temp IR Temperature of the Infrared Sensor
- Temp UV Temperature of the UV Sensor
- Pressure Gas pressure
- Rel.Humid Relative humidity
- Abs.Humid Absolute humidity
- Con.Humid Humidity concentration
- EC Sensor Electrochemical sensor (oxygen or humidity)
- O2 digi Oxygen concentration with digital O₂ Sensor



8.4. Measurement Display Refresh

The values on the display are updated automatically; no input is required.

Through the menu, values that are not relevant (e.g., because the corresponding sensor is not installed) can be hidden.

8.5. Default Settings


Unless changed, the refresh interval is preset to 10 seconds. This value can be changed in the menu. The temperature is also set to °C by default and can be displayed in °F in the menu if required. Unless changed, the menu language is set to English by default.

All changes made from the default settings are saved once the changes are confirmed by pressing the »Menu/Enter« button (blue) and will be retained after restarting the device.

8.6. Menu Navigation – Description of Menu Items

Press the »Menu/Enter« button (blue) to open the menu.

Use the "Up/Down" buttons to scroll through the individual menu items. If a menu item is activated such that the "Up/Down" buttons take over the function of value setting, this is indicated to the user by the moving indicator from the main menu item to the sub-item, and additional functions become available.



The menu includes the following items:

1. Exit:

To exit the menu, select the menu item Exit and press the »Menu/Enter« button (blue). This item is preselected when the menu is opened, so the menu can be exited immediately by pressing the »Menu/Enter« button again.



```
*MENU
EXIT
```

2. Interval

The Interval menu item allows adjustment of the measurement interval. It can be set from a few seconds up to several hours.

The value is saved and remains active after a restart.




```
*MENU
INTERVAL
60 sec
```

3. Measurement History

The History stores the last 10 timestamped measurement values for later viewing.

Procedure:

- Select the menu item »Measurement History«
- Use the »Up/Down« buttons to select the desired timestamp
- Confirm with the »Menu/Enter« button (blue)
- Scroll through the stored measurement values using the »Up/Down« buttons



```
MENU
*HISTORY
TIME: 0h1m55s AGO
```

Since the sensor unit does not have a real-time clock, the elapsed time since the measurement was recorded is shown in the format hh:mm:ss.

Note: The measurement history values are not saved after the sensor is switched off.

4. Temperature Unit

The »Temperature Unit« menu allows the display unit to be set to either °C or °F.

5. Language

The »Language« menu item allows selection of the menu language. The following languages are currently available: English, German, French, Spanish, Portuguese, and Russian (Latin alphabet).

The main display screen is not yet affected by the language settings in version v1.4. The language selection currently applies to the menu only.

6. to 12. Show/Hide Measurement Values

Menu items 6 to 12 are used to activate or deactivate optional measurement values on the main display screen, in cases where these values are not relevant for the measurement or, for example, if the corresponding sensor is not installed.



```
*MENU
SHOW CH 1/2
ACTIVE
```

13. Set Zero PPM%


The menu item »Set Zero PPM%« is used to define the zero point for channels Ch 1 to Ch 4.

14. Endpoint Correction:

The menu item »Endpoint Correction« allows calibration using a reference gas for channels Ch 1 to Ch 4.

15. Backlight:

The menu item "Backlight" allows the display backlight to be turned off. However, note that without backlighting, the display is only readable under strong external lighting. This function should therefore only be used if the display light is disruptive in the operating environment.



8.7. Menu Navigation – Overview of Menu Items

Graphical overview of the menu navigation:

